sartorius stedim
biotech

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Legal Notices | Impressum
English
Deutsch

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If you require further information on SENSOLUX® stand-alone, other than mentioned in this User Manual, please contact us.
Reproduction, translation and duplication in any form, also in part, requires written permission of the Sartorius Stedim Biotech GmbH.

Address of the technical documentation:

Sartorius Stedim Biotech GmbH
Technical editorship
August-Spindler-Strasse 11
D-37079 Goettingen, Germany
E-mail: henrik.bodenhausen@sartorius-stedim.com
Internet: www.sartorius-stedim.com

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Version: November 2009
Bestellnummer: 85034-538-96

Contents

1. SENSOLUX® stand-alone	7
1.1 Design	8
1.2 Equipment Supplied	8
1.3 Assembly of the SENSOLUX® Tray.....	9
1.4 Laying the Electrical Connection	11
2. SENSOLUX® e	14
2.1 Description of Design and Functionality	14
2.2 Electrical Connections	14
2.3 Installation and Startup	15
3. SENSOLUX® control.....	16
3.1 Description of Functionality	16
3.2 System Requirements.....	16
3.3 Installation, Equipment and Facilities.....	16
4. Entering and Running Experiments in SENSOLUX® control.....	17
4.1 Planning Example.....	17
4.2 Scanning Lot Data.....	17
4.3 Reference Method for Determining the Initial pH.....	17
4.4 Create a New Experiment.....	18
4.5 Calibration.....	21
4.6 Equilibration	25
4.7 Starting an Experiment	25
4.8 Displaying Analysis Data	26
4.8.1 Displaying in Graphs.....	26
4.8.2 Displaying in Tables	27
4.8.3 Magnification Function.....	27
4.9 Pausing an Experiment.....	28
4.10 Recalibration	28
4.11 Data Export	29
4.12 Processing Experiments.....	29
4.13 Tray Status.....	31
4.14 Deleting an Experiment.....	31
4.15 Exiting SENSOLUX® control	31
5. Appendix	32
5.1 Cleaning, Maintenance and Repairs.....	32
5.2 Warranty & Customer Service.....	32
5.3 Disposal.....	33
6. Specifications	34
6.1 SENSOLUX® tray	34
6.2 SENSOLUX® e.....	34
6.3 SENSOLUX® BCS	34
7. Ordering Information	35
7.1 Ordering Information.....	35
7.1.1 SENSOLUX® stand-alone.....	35
7.1.2 SENSOLUX® EF.....	35
7.1.3 SENSOLUX® C.....	35
7.1.4 SENSOLUX® BCS.....	35
7.1.5 Accessories.....	35

General

Notes Regarding these Instructions

This manual is part of the product. Keep it in a safe and easily accessible location.

You can request a replacement copy if you lose the manual.

This manual on SENSOLUX® stand-alone describes the product at the time of printing stated in the revision number. All data have been compiled with care, although no guaranty can be assumed for the data given.

Sartorius Stedim Biotech GmbH reserves the right to make technical changes to the design and equipment of its products or to make changes to this document without any prior notice or specific mention.

Because Sartorius Stedim Biotech GmbH 's systems are constantly undergoing further development, some components and their operation may not be described in this manual. This manual may also contain information that does not apply to your system. Please notify us if any information is missing that is important to you, if you find any errors in content or desire more extensive information on specific components.

SENSOLUX® stand-alone is just one product from the line of laboratory equipment offered by Sartorius Stedim Biotech GmbH. The rest of the product line includes shakers and shaker incubators, homogenizers, centrifuges, ultrapure water systems, air samplers, fermenters and bioreactors. For more information, please visit our website:

www.sartorius-stedim.com

On request, we would be delighted to inform you about our entire product range.

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.3118
Fax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

Symbols and Signs

The following signs and symbols are used in this manual:

These notes identify hazards that can be expected with great probability and can lead to serious injuries and possibly death.

These notes identify hazards that can occur and may lead to moderate or mild injuries.

These notes identify hazards associated with the risk of material damage.

This symbol identifies useful information and tips.

EC Declarations of Conformity

The Declarations of Conformity enclosed herein or with the SENSOLUX® stand-alone documentation provided to the customer is testimony by Sartorius Stedim Systems GmbH that its SENSOLUX® tray and SENSOLUX® e comply with the stated guidelines and directives.

Intended Use

SENSOLUX® stand-alone is an intelligent shaker tray with an integrated sensor system. It is used for monitoring and recording the pH value and the dissolved oxygen (DO) saturation during the cultivation of animal cells in shaker flasks. Any other use beyond this is considered improper.

SENSOLUX® stand-alone has been constructed in accordance with the directives and standards for operation of electrical equipment, electromagnetic compatibility and applicable safety requirements. Improper use or handling of the system, however, can result in damage and/or injury.

Follow the instructions below to ensure proper installation and operation of your SENSOLUX® stand-alone:

- This operating manual was developed according to regulations and legal guidelines currently available at the time. It has been drafted so that it can be implemented and used by appropriately trained personnel. Activities that go beyond the scope of what has been described here are to be carried out by Sartorius Stedim Biotech GmbH service personnel exclusively.
- For proper use, SENSOLUX® trays are to be operated in a CO₂ incubation shaker (e.g. Sartorius CERTOMAT® CTplus).

Due to the light-sensitive sensor patches in the flasks, you should avoid the effects of direct light on the shaker incubation room.

- Laboratory managers or persons responsible for equipment deployment are required to instruct their personnel about the hazards that can arise during the installation and operation of SENSOLUX® stand-alone, while ensuring compliance with the relevant safety rules and regulations. The personnel must possess the required qualifications for the area of application and handling of SENSOLUX® stand-alone.
- Strictly follow and observe the installation and operating instructions. That way you can avoid improper handling and maintain all guarantee claims in full. Please contact us if you have any questions.
- If these safety instructions are not followed, this can cause damage to your SENSOLUX® system or to other equipment and material or even personal injury.
- Every time before using the equipment, check that all components are in perfect working order.

Never use damaged components.

- Observe the max. allowable load for the shaker and ensure that there is enough space available for the components.
- Only operate the SENSOLUX® tray within the allowable operating tolerances (see Chapter 6.1 SENSOLUX® tray).

Guide to these Instructions

The operating instructions manual is set up as described below. This overview will help you to directly find the information that is important to you and that satisfies your information needs.

Chapter 1

describes the SENSOLUX® stand-alone design and functionality as well as how to install the SENSOLUX® tray in an incubation shaker.

Chapter 2

describes the SENSOLUX® e data collector unit and its function as well as how to connect the SENSOLUX® tray, computer system and barcode scanner to the data collector unit.

Chapter 3

describes the SENSOLUX® control software for planning and conducting experiments, monitoring and visualizing measured data as well as for installing the system on a computer system (PC/laptop).

Chapter 4

describes how to program and conduct experiments with the SENSOLUX® control software as well as options for exporting and/or further processing the recorded data.

Chapter 5

lists cleaning, maintenance instructions and information about the warranty and customer service.

Chapter 6

provides an overview of the specifications for the SENSOLUX® stand-alone.

Chapter 7

contains ordering information. This is where you will find special information about optional equipment and the accessories program offered with the SENSOLUX® stand-alone.

1. SENSOLUX® stand-alone


SENSOLUX® Technology

SENSOLUX® stand-alone is an intelligent shaker tray featuring an integrated sensor system. It is intended for monitoring pH values and dissolved oxygen (DO) saturation during the cultivation of mammalian cells in shake flasks.

The non-invasive measuring principle for the two analytes is based on an optical procedure. Each single SENSOLUX® Erlenmeyer Flask (EF) contains pre-calibrated sensor patches for pH and DO. The SENSOLUX® shaker tray features 9 independent measurement points (for pH and DO, respectively) for reading out these sensor patches.

The sensor patches at the bottom of the flasks contain fluorescent dyes, which can be excited to generate fluorescence by light of a given wavelength shining through the flask wall. The properties of the longer-wavelength fluorescence light emitted by the sensor patches are respectively influenced by the parameters pH or dissolved oxygen saturation of the culture medium in contact with the sensor patches. The SENSOLUX® shaker tray is equipped with optoelectronics that run both the fluorescence excitation and the capture and analysis of the fluorescence light.

The optoelectronics transmit the measured values for DO and pH to the data collector unit (SENSOLUX® e) and the software (SENSOLUX® control) for processing. Used in combination with the new ready-to-use disposable SENSOLUX® Erlenmeyer flasks, the SENSOLUX® shaker tray facilitates the easy, safe and highly informative online measurement of these crucial process parameters in incubation shakers during cultivation.


Features & Benefits of SENSOLUX® Technology

- Non-invasive measurement
- Online detection
- Fast
- Reproducible
- Ready-to-use

1.1 Design

SENSOLUX® stand-alone consists of a shaker tray equipped with an optical sensor system. It is applied for monitoring the pH value and the dissolved oxygen (DO) saturation during the cultivation of animal cells. Used in combination with the new single-use SENSOLUX® Erlenmeyer flasks, it enables non-invasive measurement in incubation shakers. A total of nine measurement points for detecting pH and DO are installed in the shaker tray.

Schematic diagram of SENSOLUX® stand-alone and its components

1.2 Equipment Supplied

SENSOLUX® stand-alone consists of the following components:

- SENSOLUX® tray – Shaker tray with integrated optical sensors (9 measurement points) for monitoring pH value and dissolved oxygen (DO) saturation in shake flasks,
- Adapter plate for installing SENSOLUX® trays on the shaker platform,
- SENSOLUX® e – Data collector unit,
- SENSOLUX® control – Software (on CD/DVD) for planning and conducting experiments as well as for monitoring and visualizing the measured data,
- Full set of cabling required for installing SENSOLUX® stand-alone:
 - Connecting cable SENSOLUX® D1
(Cable for connecting SENSOLUX® tray and incubator wall),
 - Connecting cable SENSOLUX® D2
(Cable for connecting SENSOLUX® D1 and SENSOLUX® e),
 - Connecting cable SENSOLUX® D3
(Cable for connecting SENSOLUX® e and PC/laptop),
 - Power supply with power supply cable (230V/115V) for SENSOLUX® e
- Optional:
 - SENSOLUX® BCS – Barcode scanner for automatic and safe transfer of the calibration data from pH and DO sensor patches.

1.3 Assembly of the SENSOLUX® Tray

- Remove tray and shaker platform from your incubation shaker.
- Attach the adapter plate supplied by Sartorius to the mechanical shaker parts of your incubation shaker.
 - For assembling the adapter plate, use the original fixing screws of the shaker platform.

When attaching the adapter plate to the mechanical shaker parts, pay attention to the proper orientation of the adapter plate.

Proper orientation means that the hole pattern of the adapter plate matches up with the hole pattern of the mechanical shaker parts.

- Next, attach the shaker clamps (not included with the standard equipment supplied) to the SENSOLUX® tray.

The following types of shaker clamps can be attached to the SENSOLUX® tray:

Type	Description
SENSOLUX® C 125	SENSOLUX® EF 125 ml stainless steel clamps for single-use Erlenmeyer flasks
SENSOLUX® C 250	SENSOLUX® EF 250 ml stainless steel clamps for single-use Erlenmeyer flasks
SENSOLUX® C 500	SENSOLUX® EF 500 ml stainless steel clamps for single-use Erlenmeyer flasks
SENSOLUX® C 1000	SENSOLUX® EF 1000 ml stainless steel clamps for single-use Erlenmeyer flasks


We recommend that you attach the shaker clamps to the SENSOLUX® tray before you insert the tray in the incubation shaker. This method allows for easier handling.


During assembly, make sure that the shaker clamps are evenly distributed around the SENSOLUX® tray. This method will best keep the tray balanced during operation later.

- Attach the shaker clamps so that the recessed spaces in the base plate are pointing to the back (see the following drawing).


Important Note:
Never cover the sensor system with shaker clamps!

- Hand-tighten the fixing screws evenly and firmly (see following drawing).

Proper assembly

Incorrect assembly

- Attach the loaded SENSOLUX® tray to the adapter plate attached to the incubation shaker.
- For secure assembly, use the 8 hexagon socket screws (included with the equipment supplied).


1.4 Laying the Electrical Connection


On incubation shakers factory-equipped with cable entry, first draw the black ribbon cable SENSOLUX® D2 into the incubation chamber so that the right-angle plug for connection to the SENSOLUX® tray is located in the incubation chamber. The connection to SENSOLUX® e is located outside of the incubation chamber.

- For fastening the right-angle plug in the incubation chamber, place it in the intended mounting angle.
- Connect the white ribbon cable SENSOLUX® D1 to the right-angle plug and secure it with the plug screws.
- Connect the other end of the white ribbon cable SENSOLUX® D1 to Port 1 of the SENSOLUX® tray.
- Hold the mounting angle up against the intended mounting position (Do not affix yet!).

- Turn the shaker motor to ensure that the ribbon cable has enough room in every position during orbital movements of the SENSOLUX® tray.

- Remove the protective film of the adhesive pad on the mounting angle and glue this to the back panel of the incubation chamber.
- When gluing in the mounting angle, pay attention:
 - That adhesive surfaces are clean and free of grease,
 - To the lateral and height alignment of the mounting angle (see following drawing).

An incorrectly mounted right-angle plug can cause excess wear to the white ribbon cable SENSOLUX® D1.

- Lay the black ribbon cable SENSOLUX® D2 in the incubation chamber of your shaker.

Sartorius recommends laying the cable in the ceiling area of the incubation chamber.

- For fastening the ribbon cable, use the cable holders included in the equipment supplied by Sartorius.
Affix the cable holders to the appropriate sites in the incubation chamber.
- Make sure that the adhesive surfaces are clean and free of grease.

When feeding the ribbon cable through the door:

- In the area of the door feed-through, pay attention to:
 - Door contacts,
 - Hinges,
 - Sharp edges.

These may damage the ribbon cable.

2. SENSOLUX® e

2.1 Description of Design and Functionality

The data collector unit (SENSOLUX® e) serves as the interface between the SENSOLUX® control software and the SENSOLUX® tray and enables them to communicate with each other. Moreover, the data measured by the SENSOLUX® tray are stored there. Even in the event the computer (PC/laptop) crashes, the stored experiments and data remain secured on the data collector unit (SENSOLUX® e).

2.2 Electrical Connections

Data collector unit on the SENSOLUX® e

1. ON/OFF switch
2. USB port (for optional barcode scanner)
3. Network port (currently nonfunctional)
4. Port for ribbon cable SENSOLUX® D2
5. Sartorius Service interface
6. Port for SENSOLUX® D3 (connection to PC/laptop)
7. Port for power supply cable

2.3 Installation and Startup

Do not connect the data collector unit SENSOLUX® e to your PC/laptop until after you have installed the SENSOLUX® control software (see Chapter 3 SENSOLUX® control).

Do not switch the data collector unit on using the ON/OFF switch until you have connected all cable.

To start up the data collector unit SENSOLUX® e, proceed as follows:

- Connect the ribbon cable SENSOLUX® D2 (connection to SENSOLUX® tray) to the data collector unit.
 - Connect the connecting cable SENSOLUX® D3 (connection to PC/laptop) to the data collector unit.
 - Connect the other end of the cable to a free USB port on your PC/laptop.
 - Connect the optionally available barcode scanner SENSOLUX® BCS to the data collector unit.
 - Connect the power supply cable to the data collector unit.
 - Switch the power supply on. As per your order, the power supply includes a country-specific connecting cable.
 - Switch on the data collector unit. The ON/OFF switch is located on the front panel of the data collector unit.
-

After switching on the data collector unit, it will run an initialization routine. The initialization procedure takes approx. 30 seconds. Afterwards, the data collector unit is operational.

3. SENSOLUX® control

3.1 Description of Functionality

SENSOLUX® control software is used for planning and conducting experiments as well as for monitoring and visualizing the measured data. Moreover, the recorded data can be exported in an XML format (compatible with Excel) and are available for further processing.

3.2 System Requirements

To install and run SENSOLUX® control on your PC/laptop, your system must meet the following requirements:

Processor type	AMD/Intel 32 bit processor, min. 2 GHz
Memory	min. 512 MB RAM (1 GB are recommended)
Free hard drive space	min. 1 GB <u>are</u> recommended
Graphic card	min. resolution 1024 x 768 or higher
Operating system	Microsoft Windows XP with Service Pack 2 or higher; Microsoft Vista; Microsoft Windows 7
Additional software components	Microsoft .net Framework 3.5 SP1

3.3 Installation, Equipment and Facilities

- To install SENSOLUX® control, place the installation CD/DVD in the CD/DVD drive.
- The installation procedure starts automatically.
- During installation, a check is run to make sure that all the software components are installed on your PC/laptop that are necessary for running SENSOLUX® control. Under certain circumstances, administrative rights may be necessary to install additional components.
- Follow the further installation instructions.
- After successful installation, switch on the data collector unit (SENSOLUX® e).

Note that the initialization time of the data collector unit is approx. 30 seconds.

- Connect the data collector unit (SENSOLUX® e) to a free USB port on your PC/laptop.
- Your PC/laptop will automatically identify and install the data collector unit (SENSOLUX® e).
- Start SENSOLUX® control by double-clicking on the desktop icon.

4. Entering and Running Experiments in SENSOLUX® control

4.1 Planning Example

The following data are chosen for the example experiment:

Experiment	Test
Measurement interval	1 minute
Description (optional)	Sart test
Number of flasks	5

4.2 Scanning Lot Data

The lot data can be scanned at any time using the barcode scanner as long as the data collector unit (SENSOLUX® e) is operational. This is independent of SENSOLUX® control.

- Use the barcode scanner to scan the barcode on the sticky label affixed to the cardboard box with the flasks. The barcode contains the calibration data on the box containing the flasks.

Lot data label (affixed to the box containing the flasks)

You can pull-off the lot data label from the box (peel-off label) and save it for internal documentation.


4.3 Reference Method for Determining the Initial pH

The following method is suited for determining the pH of the medium before starting cultivation:

- Working under a clean bench, fill approx. 1 ml of medium into a sealable sample vessel (e.g. test tube),
- Place the unsealed sample vessel into the CO₂ incubator until the temperature and CO₂ equilibrium are reached (approx. 20-30 minutes),
- Before removing the sample vessel, seal it,
- Shortly before immersing the pH electrode in the sample vessel, open the latter and measure the pH of the medium (e.g. using the Basic pH meter PB-11 supplied by Sartorius Mechatronics).


4.4 Create a New Experiment

- Start SENSOLUX® control by double-clicking on the desktop icon.

Main overview with all stored experiments and their respective status.

- To create a new experiment, click on "Experiment" --> "New" in the header.


Input screen for creating a new experiment.

- Enter the name of the experiment in the "Experiment" field.
- Open the "Measurement interval" field with a click.
- Select one of the pre-set intervals.
- In the "Description" field, you have the option to enter a brief description of the experiment.
- Deactivate all flasks you do not need with a click.
 - They are displayed by colored shading (inactivated).

Ensure even distribution on the SENSOLUX® tray.
This method will best keep the tray balanced during operation later.


Optionally, you can label the flasks you need by giving them a name/or affixing a label to them.

- To do so, click in the white labeling area of the activated flasks and label/name this accordingly.

When labeling the flasks, pay attention that each flask has a different name.
If the names of the flasks are the same, it will not be possible to create an experiment.

- Click on "Ok", to confirm your entries.

The main view shows all previously stored experiments and the newly created experiment named "Test".

- Fill the flasks with your medium for equilibration.

Please consider:

The flasks are individually packaged with light protection. The sensor patches are very sensitive to light, therefore the flasks should only be exposed to any light sources as briefly as possible.

If the desired end volume should be x ml, please pay attention to the volume of the cell suspension to be entered later.

End volume x = charged medium + cell suspension

4.5 Calibration

Before starting an experiment, the calibration data of the flasks and the "Initial pH" value have to be scanned or entered in the calibration menu.

- Press the "Calibration" button for the experiment to be calibrated.

Experiment

Please choose your flask below to set calibration values for all of them at once.

Flask 1

Initial pH [-] Initial pO2 [%] Temp. [°C]

pH [-] pmin pmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Flask 3

Initial pH [-] Initial pO2 [%] Temp. [°C]

pH [-] pmin pmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Flask 5

Initial pH [-] Initial pO2 [%] Temp. [°C]

pH [-] pmin pmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Flask 7

Initial pH [-] Initial pO2 [%] Temp. [°C]

Cancel Ok

Menu for inputting calibration data and "Initial pH".

- Select the flasks for which you would like to enter a mutual basic calibration by inactivating the flasks not needed with a click.

Experiment

Please choose your flask below to set calibration values for all of them at once.

Flask 3

Initial pH [-] Initial pO2 [%] Temp. [°C]

pH [-] rmin rmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Flask 5

Initial pH [-] Initial pO2 [%] Temp. [°C]

pH [-] rmin rmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Flask 7

Initial pH [-] Initial pO2 [%] Temp. [°C]


pH [-] rmin rmax dpH pH0

pO2 [%] Phi (100%) Phi (0%)

Cancel X Ok ✓

Inactivated flasks are displayed with colored shading.

- Take a pre-filled flask after you have already scanned in its lot data (see Chapter 4.2), and scan the barcode affixed to the flask. These data will be compared with the lot data of the already stored lots, thereby enabling an unequivocal assignment of the flasks to the respective lot and its calibration data. The lot data will be transferred to the input screens of the selected flasks of the calibration menu.

Lot No. Label (affixed to the flask)

- Place the flasks in the shaker clamps on the SENSOLUX® tray according to the experiment as programmed in the software. The bottom sides of the flasks have positioning feet. These feet must lock into the holes on the shaker clamps; the correct orientation of the flasks is noted by tactile locking-in.

The mechanism ensures the exact positioning of the sensor patches and thereby precise measurement.

- Enter initial pH.
- Confirm entry with "Tab key".

- As an alternative, whenever there is no barcode scanner available for identifying the flasks, you can also manually enter the data contained on the box label (for lot data, see Chapter 4.2) in the input screen,
- When manually entering data, use the "Tab key" to move around from field to field within the input screen.
Confirm entry of last value by pressing "Tab key".

The figure shows the flasks calibrated in a lot displayed in yellow.

- Activate the remaining flasks with a click to scan or enter their lot data.
- Take a pre-filled flask after you have already scanned in its lot data (see Chapter 4.2), and scan the barcode affixed to the flask.
These data will be compared with the lot data of the already stored lots, thereby enabling an unequivocal assignment of the flasks to the respective lot and its calibration data. The lot data will be transferred to the input screens of the selected flasks of the calibration menu.

Lot No. Label (affixed to the flask)

- Place the flasks in the shaker clamps on the SENSOLUX® tray according to the experiment as programmed in the software. The bottom sides of the flasks have positioning feet. These feet must lock into the holes on the shaker clamps; the correct orientation of the flasks is noted by tactile locking-in.

The mechanism ensures the exact positioning of the sensor patches and thereby precise measurement.

- Enter initial pH.
- Confirm entry with "Tab key".
- As an alternative, whenever there is no barcode scanner available for identifying the flasks, you can also manually enter the data contained on the box label (for lot data, see Chapter 4.2) in the input screen,
- When manually entering data, use the "Tab key" to move around from field to field within the input screen.
Confirm entry of last value by pressing "Tab key".

The figure shows the selected and calibrated flasks displayed in yellow.

- Confirm calibration by pressing "Ok" button.

4.6 Equilibration

The SENSOLUX® control software presets the equilibration time at 3 hours. You cannot change it.

- Toggle, if necessary, to the main view.
- Press the "Equilibration" button for the experiment to be equilibrated.
- Next, press the "OK" button.

The experiment cannot be started before equilibration time (3 hours) has elapsed.

During the equilibration time, you can create, process or start new or existing experiments (depending on the status), or pause or delete them.
After the equilibration time has elapsed, you can inoculate the individual flasks.

Keep the period between sampling, inoculation and returning the individual flasks to the incubation shaker as short as possible.

4.7 Starting an Experiment

After adjusting the inoculated flasks, the inoculated medium must return to a state of equilibrium (temperature and CO₂).

- Toggle, if necessary, to the main overview.
- Click on the "Start" button for the experiment to be started.
- Next, press the "OK" button.


4.8 Displaying Analysis Data

You can display the previously obtained measured values from all stored experiments.

- Toggle, if necessary, to the main overview.
- With a click, select the experiment for which you would like to display the existing measured values.
- The selected experiment is highlighted in color.

4.8.1 Displaying in Graphs

To display the measured values as graphic output, press the "F5" function key or click on "View" --> "Graphic".

The figure shows the curve of the measured values for the experiment entitled "Test" plotted over a period of 1 hour for the individual flasks.

- Use the "Time scale" slide bar to set a suitable value for the respective experiment.

4.8.2 Displaying in Tables

- To display the measured values in table format, press the "F6" function key or click on "View" --> "Data".


The figure shows the tables of the measured values from the experiment "Test" for the individual flasks.

4.8.3 Magnification Function


As needed, for the values of an individual flask, you can generate a magnified display of the curve in the graph view (F5) and/or of the table of values in the data view (F6).

- Move the mouse pointer to the upper right-hand corner (see picture below) of the flask, its curve, and/or table of values of which you want to magnify the display.
- The magnification icon will appear.
- Click on the icon to toggle to the magnified display.
- To toggle back to the normal view, press F5 (graph view) or F6 (table view) or move the mouse pointer to the upper right-hand corner of the magnified view.
- The minimization icon will appear.
- Click on the icon to switch to the normal view.

Example: Graph view

Magnification icon

Minimization symbol

4.9 Pausing an Experiment

- Toggle, if necessary, to the main overview.
- Click on the "Pause" button for the experiment to be paused.
- Next, press the "OK" button.

4.10 Recalibration

If, during an already started experiment, you need to change experiment data on individual or all flasks, you can do so using the "Recalibration" function.

- Toggle, if necessary, to the main overview.
- Click on the "Pause" button for the experiment to be recalibrated.
- Next, press the "OK" button.
- Click on the "Recalibration" button for the experiment to be recalibrated.

Figure shows the input menu for the recalibration.

- Enter new values in the value fields for the flasks whose values have changed.
- If the values of the individual flasks deviate from one another, remove the checkmark in the field "Update" for the respective flask.
- Click "All", if you want to transfer the changed value to all flasks of the experiment.
- Exit the entry by pressing "OK" button.
- Next, after you have properly placed all flasks back on the tray, press the "OK" button.

The experiment starts with the updated values.

4.11 Data Export

The "Excel XML format" is available for exporting data.
You can process and analyze exported data in all programs compatible with this format.

- To export data, select the experiment to be exported with a click in the Main overview; in the header, select "Experiment" --> "Export".
- Assign a file name and a memory path and confirm with "Save".

4.12 Processing Experiments

You can change the basic data of an experiment with the "Processing" menu item.

- Experiment
- Measurement interval
- Description
- Number of flasks (only in the "Planned" or "Pause" status)

Status "Planned": Number of flasks can be increased or decreased.
Status "Pause": Number of flasks can only be decreased.

To change basic data in the experiment status "Pause" proceed as follows:

- Toggle, if necessary, to the main view.
- Select the experiment to be processed (Status "Pause") with a click.
 - The selected experiment is highlighted in blue.
- Click in the header on "Experiment" --> "Processing".

Basic data of the experiment to be processed

- Change the basic data as required.

Important Note:
It is not possible to activate deactivated flasks!

- Click on the flask(s) to be deactivated.


Overview of the processed basic data indicating deactivated flasks

- Exit entries by pressing "OK" button.

4.13 Tray Status

The "Tray Status" menu item provides information about the connection between the data collector unit and the shaker tray. The "Tray Status" additionally gives you the option for reinstating the connection after interruption of the connection.

- To display the tray status, click on "View" --> "Tray Status" in the header.

- To exit this view, click on "Back".

4.14 Deleting an Experiment

Important Note: Deleted experiments are deleted irrevocably.

- To delete an experiment, click on it in the main overview, select "Experiment" --> "Delete" in the header and then confirm with "Ok".

4.15 Exiting SENSOLUX® control

- To exit SENSOLUX® control, click on "Experiment" --> "Exit program" in the header.

Data from experiments with status "Active" before you exited SENSOLUX® control will continue to be recorded and saved in the SENSOLUX® e data collector unit.

The next time SENSOLUX® control is started, the data will be available for visualization and/or further processing.

5. Appendix

5.1 Cleaning, Maintenance and Repairs

You can clean dirty or contaminated surfaces on the SENSOLUX® tray with conventional household cleaning agents or alcohol.

Avoid aggressive cleaning agents.
Agents containing chlorine, for example, can cause corrosion

When cleaning, pay attention that you do not damage the optical window on the SENSOLUX® tray.

All components of SENSOLUX® stand-alone are maintenance-free.

Repairs on the SENSOLUX® tray (e.g. the replacing of defect sensor system elements or other electrical components on the SENSOLUX® tray) are only allowed to be performed by qualified and authorized service technicians.

5.2 Warranty & Customer Service

Unless agreed otherwise in writing, all equipment from Sartorius Stedim Biotech GmbH is covered by a warranty in accordance with the General Terms and Conditions of Purchase as per the date of delivery.

- The warranty applies to defects and defective operation caused by faulty construction, manufacturing and material faults, but not for disturbances that are attributable to improper operation and improper handling and use. The warranty also does not apply to parts that are defective due to normal wear and tear or to consumables.
- The warranty is invalidated whenever the user or unauthorized third parties make technical changes to SENSOLUX® stand-alone and/or use equipment, components and accessories that have not been approved by Sartorius Stedim Biotech GmbH for use in conjunction with SENSOLUX® stand-alone.

Your Sartorius Stedim Biotech GmbH service dealership or local authorized service technicians are happy to repair any defects for you.

In the event of a defect, please notify your Sartorius Stedim Biotech GmbH dealer or contact directly:

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
info@sartorius-stedim.com
www.sartorius-stedim.com

- When performing repairs, only replacements may be used that have been approved by Sartorius Stedim Biotech GmbH for use in conjunction with SENSOLUX® stand-alone.
- You can return defective components to Sartorius Stedim Biotech GmbH.

- Sartorius will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal.
Please refer to the package leaflet | the included CD ROM or to our web site (www.sartorius-stedim.com) for more detailed information and addresses for repair service or disposal of your device.

Returned components must be clean and in hygienically flawless condition and packed carefully. Parts contaminated by media or media ingredients must be disinfected, decontaminated and/or sterilized in accordance with the safety guidelines applicable to the area application, e.g. those concerning chemical or biological safety.

The sender must prove compliance with the safety guidelines. The enclosed Decontamination Declaration must be completed and signed and included with the components. Without proof of harmlessness or a description of the measures taken, Sartorius Stedim Biotech GmbH cannot repair the components.

Transport damage as well as subsequent cleaning and disinfection of the parts – if necessary – shall be charged to sender.

5.3 Disposal

Attention must be paid to the applicable laws, rules and regulations when disposing of the contaminated components.

If the packaging is no longer needed, it can be disposed by taking it to a local waste disposal site. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The tray – including accessories – should not be disposed of as regular household waste.

EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and several other countries, Sartorius Stedim Biotech GmbH itself assumes responsibility for the return and conformant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH
Biotechnology Service Center
August-Spindler-Strasse 11
37079 Goettingen, Germany
WEEE Reg No. DE 89907997

In countries that are not members of the European Economic Area (EEA) or where no Sartorius Stedim Biotech subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and | or scrapping of the components, any batteries – if present – should be removed and disposed of in local collection boxes.

6. Specifications

6.1 SENSOLUX® tray

Dimensions (W + D + H)	406 x 420 x 62.5 mm
Weight	10 kg
Material	Stainless steel/Aluminum/Pyrex
Interface ports/connections	SENSOLUX® Bus (DSUB-9, female)
Number of measurement points	9
Max. revolutions per minute	250 rpm, during rotational movement Ø 50 mm
Max. load	8 kg
Ambient temperature	+5 °C ... +45 °C
Humidity	0 % ... 99% (non-condensing)

6.2 SENSOLUX® e

Dimensions	160x50x184 mm
Weight	620g
Interface ports/connections	USB RS232 (DSUB-9, female) SENSOLUX® Bus (DSUB-9, male) Ethernet (RJ45)
Supply voltage	24 VDC, 2.2 A
Ambient temperature	+5 °C ... +45 °C
Humidity	0 % ... 95 % (non-condensing)

6.3 SENSOLUX® BCS

Interface ports/connections	USB
Supply voltage	5 VDC (per USB)
Ambient temperature	0 °C ... +50 °C
Humidity	0 % ... 95 % (non-condensing)

7. Ordering Information

7.1 Ordering Information

7.1.1 SENSOLUX® stand-alone

	Product Description	Art. No.
SENSOLUX® stand-alone	The SENSOLUX® stand-alone – shaker tray with a sensor system (9 measurement points) for optical pH and DO determination in shake flasks, incl. set of cables with male connectors, data collector unit (SENSOLUX® e) and software (SENSOLUX® control).	DCS09

7.1.2 SENSOLUX® EF

	Product Description	Art. No.
SENSOLUX® EF 125	Single-use 125 ml Erlenmeyer flasks with a vented cap, pH and DO sensors, 9 pcs/case	DCS--F1
SENSOLUX® EF 250	Single use 250 ml Erlenmeyer flasks with a vented cap, pH and DO sensors, 9 pcs/case	DCS--F2
SENSOLUX® EF 500	Single-use 500 ml Erlenmeyer flasks with a vented cap, pH and DO sensors, 9 pcs/case	DCS--F3
SENSOLUX® EF 1000	Single use 1000 ml Erlenmeyer flasks with a vented cap, pH and DO sensors, 9 pcs/case	DCS--F4

7.1.3 SENSOLUX® C

	Product Description	Art. No.
SENSOLUX® C 125	SENSOLUX® EF 125 stainless steel clamps, 9 pcs/case	DCS--C1
SENSOLUX® C 250	SENSOLUX® EF 250 stainless steel clamps, 9 pcs/case	DCS--C2
SENSOLUX® C 500	SENSOLUX® EF 500, stainless steel clamps, 9 pcs/case	DCS--C3
SENSOLUX® C 1000	SENSOLUX® EF 1000 stainless steel clamps, 9 pcs/case	DCS--C4

7.1.4 SENSOLUX® BCS

	Product Description	Art. No.
SENSOLUX® BCS	Barcode scanner for the automatic transfer of calibration data	DCS--B1

7.1.5 Accessories

	Product Description	Art. No.
SENSOLUX® D1	Connecting cable between SENSOLUX® tray and incubator wall	DCS09D1
SENSOLUX® D2	Connecting cable between SENSOLUX® D1 and SENSOLUX® e	DCS09D2
SENSOLUX® D3	Connecting cable between SENSOLUX® e and PC/Laptop	DCS09D3
SENSOLUX® D4	Power supply cable and power supply (230V/115V) for SENSOLUX® e	DCS09D4

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For further contacts, visit www.sartorius-stedim.com

For further contacts, visit www.sartorius-stedim.com

Germany
Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

France
Sartorius Stedim Biotech S.A.
ZI Les Paluds
Avenue de Jouques – BP 1051
13781 Aubagne Cedex
Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI Les Paluds
Avenue de Jouques – CS 71058
13781 Aubagne Cedex
Phone +33.442.845600
Fax +33.442.846545

Austria
Sartorius Stedim Austria GmbH
Franzosengraben 12
A-1030 Vienna
Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium
Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde
Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark
Sartorius Stedim Nordic A/S
Hoerskaetten 6D, 1.
DK-2630 Taastrup
Phone +45.7023.4400
Fax +45.4630.4030

Italy
Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)
Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands
Sartorius Stedim Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein
Phone +31.30.6025080
Fax +31.30.6025099

Spain
Sartorius Stedim Spain SA
C/Isabel Colbrand 10-12,
Planta 4, Oficina 121
Poligono Industrial de Fuencarral
28050 Madrid
Phone +34.91.3586102
Fax +34.91.3588804

Switzerland
Sartorius Stedim Switzerland GmbH
Lerzenstrasse 21
8953 Dietikon
Phone +41.44.741.05.00
Fax +41.44.741.05.09

U.K.
Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ
Phone +44.1372.737159
Fax +44.1372.726171

USA
Sartorius Stedim North America Inc.
5 Orville Drive
Bohemia, NY 11716
Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Sartorius Stedim SUS Inc.
1910 Mark Court
Concord, CA 94520
Phone +1.925.689.6650
Toll Free +1.800.914.6644
Fax +1.925.689.6988

Sartorius Stedim Systems Inc.
201 South Ingram Mill Road
Springfield, MO 65802
Phone +1.417.873.9636
Fax +1.417.873.9275

Argentina
Sartorius Argentina S.A.
Int. A. Avalos 4251
B1605ECS Munro
Buenos Aires
Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil
Sartorius do Brasil Ltda
Av. Dom Pedro I, 241
Bairro Vila Pires
Santo André
São Paulo
Cep 09110-001
Phone +55.11.4451.6226
Fax +55.11.4451.4369

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China
Sartorius Stedim Beijing
Representative Office
No. 33, Yu'an Road,
Airport Industrial Zone B, Shunyi District
Beijing 101300
Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Shanghai
Representative Office
Room 618, Tower 1, German Centre,
Shanghai, P.R.C., 201203
Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180
Phone +86.20.8351.7921
Fax +86.20.8351.7931

India
Sartorius Stedim India Pvt. Ltd.
10, 6th Main, 3rd Phase Peenya
KIADB Industrial Area
Bangalore - 560 058
Phone +91.80.2839.1963|0461
Fax +91.80.2839.8262

Japan
Sartorius Stedim Japan K.K.
KY Building, 8-11
Kita Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001
Phone +81.3.3740.5407
Fax +81.3.3740.5406

Malaysia
Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur
Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore
Sartorius Stedim Singapore Pte. Ltd.
10, Science Park Road, The Alpha
#02-25, Singapore Science Park 2
Singapore 117684
Phone +65.6872.3966
Fax +65.6778.2494

Australia
Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175
Phone +61.3.8762.1800
Fax +61.3.8762.1828